



Oakridge Co-op Mall Professional Wing
 Suite 49, 2580 Southland Drive SW
 Calgary, AB T2V 4J8
OakridgeDentureCentre@gmail.com
www.OakridgeDenture.com
 403-251-1522

Charles Gulley, DD, F.C.A.D.
 Denturist

Cleaners Tested on Partial Denture Frameworks

If you have a partial denture which has a Titanium framework rather than a Chrome-Cobalt framework, please be advised that some commercially available cleaners will react with the Titanium.

Below is a list of typical denture cleaners and whether they are safe or not to use with your Titanium partial. In general, it is not recommended to use effervescent products with Titanium frameworks.

Cleaner Name	No Reaction	Reaction
Nu-Dent™	✓	
Nitradine™ Denture Tabs	✓	
Renew™	✓	
Dentu Creme™	✓	
Protech™	✓	
Equate™	✓	
Stain-Away Plus™	✓	
Dent-Appeal™	✓	
Nova-Dent™	✓	
Polident Fresh Cleanse™	✓	
Equate Double Action™		✗ <i>Yellow Tarnish</i>
Efferdent™		✗ <i>Yellow Tarnish</i>
Kleenite™		✗ <i>Dark Blue Reaction</i>